Ultradent[®]
Clone[™] Bite

Description/Indications:

Ultradent[®] Clone[™] Bite is an orange colored addition vulcanizing silicone for all bite registration requirements with a very high final hardness and is non-slumping. It is also a fast setting material perfect for preliminary impressions for crown and bridge procedures and temporaries.

Technical Data:

- Working time is 44 seconds.
- Setting time in mouth is 60 seconds.
- Material is completely set in 90 seconds from mix.
- Final hardness is 95 Shore A.

Procedure:

1. Load cartridge in dispensing gun.
2. Immediately prior to taking an impression, remove the cartridge cap and discard.
3. Gently squeeze dispenser to bleed enough material from the cartridge to ensure an even flow of the base and catalyst, (about 1/4 of an inch).
4. Attach the mixing tip with or without intra-oral tip. Intra-oral tip is available if the operator prefers a change in angulation of the material being dispensed.
5. Squeeze the dispenser to extrude the product.
6. If taking a bite registration, ensure the patient is brought into maximum intercuspation during the initial working time of 30 seconds.
7. Once maximum intercuspation has been achieved, have the patient remain in this position for 50 seconds to ensure a total set of the material.
8. Remove the material from the patient's mouth. Due to a linear dimensional change of <0.1% the practitioner can rest assured the material is highly accurate.
9. After use, store the cartridges with used mixing tips attached to protect the remaining materials from contamination. The material can be disinfected immediately.

Precautions:

1. Due to rapid polymerization, the patient may experience a physiologically harmless temperature increase.
2. Never attempt to extrude the material through a cured mixing canula.
3. Always disinfect cartridges between patients.
4. Store at room temperature.
5. Latex or Sulfur may affect the set of material.
6. Clone Bite may stain clothing.

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Ultradent cartridges have an expiration date and lot number printed on the label.

The lot number is a five digit number used for manufacturing purposes.

The expiration date will read YYYY-MM. Store at room temperature. Keep out of reach of children.


Shelf life 36 months.

U.S. Federal law restricts this device to sale by or on the order of a state licensed dental professional.


For immediate reorder and/or complete descriptions of Ultradent's product line, refer to Ultradent's catalog or call Toll Free **1-800-552-5512**. Outside U.S. call (801) 572-4200 or visit our website at www.ultradent.com.

Ultradent® Clone™ Bite

Flammability

MSDS Required and Included.

Hazard Rating
 4 = Severe
 3 = Serious
 2 = Moderate
 1 = Slight
 0 = Minimal

Store at Room Temperature

ULTRADENT INTERNATIONAL, INC. MATERIAL SAFETY DATA SHEET

Controlled Product No Yes

SECTION I - IDENTIFICATION OF THE SUBSTANCE/PREPARATION AND OF THE COMPANY/UNDERTAKING			
Material Name: Ultradent® CloneBite		Manufactured by: Ultradent Products, Inc. 505 West 10200 South South Jordan, Utah 84095 USA	
Material Use: Bite Registration Material		Emergency Contact: 1-800-552-5512 or 1-801-572-4200 (Outside USA)	
Date: January 5, 2005		EC Representative:	
Prepared by: Peter Allred			
SECTION II - COMPOSITION/INFORMATION ON INGREDIENTS			
Active Ingredients	Approximate Concentration %	LD ₅₀ /CL ₅₀ Specify Species and Route	C.A.S. N.A. or U.N. Numbers
N/A			
SECTION III - HAZARDS IDENTIFICATION			
Hazardous Properties: N/A Effects and Symptoms: N/A			
SECTION IV - FIRST AID MEASURES			
Inhalation	Supply fresh air. Consult doctor with any complaints.	Eye Contact	Raise opened eye for several minutes with running water. Consult doctor if symptoms persist.
Skin Contact	Generally not an irritant.	Ingestion	Induce vomiting if large amounts are swallowed. Seek medical attention.
SECTION V - FIRE-FIGHTING MEASURES			
Flammability: No <input type="checkbox"/> Yes <input type="checkbox"/>	if Not Usable When Conditions: CO ₂ powder, or water spray.	Explosion	Protective Equipment for Firefighters: Standard equipment.
SECTION VI - ACCIDENTAL RELEASE MEASURES		SECTION VII - HANDLING AND STORAGE	
Personal & Environmental Precautions	Avoid contact with skin and eyes.	Handling Precautions	Always wash hands after handling.
Methods for Cleaning	Absorb with liquid-binding materials.	Storage Requirements	Keep containers tightly sealed and store in a cool location.
SECTION VIII - EXPOSURE CONTROLS/PERSONAL PROTECTION			
Respiratory Protection	Not required.	Eye Protection	Protective eyewear.
Hand Protection	Latex gloves.	Skin Protection	Lab coat.
SECTION IX - PHYSICAL & CHEMICAL PROPERTIES			
Appearance	Orange, viscous liquid.	Stability in Water	None.
Date	N/A	Other Data	None.
SECTION X - STABILITY AND REACTIVITY		SECTION XI - TOXICOLOGICAL INFORMATION	
Chemical Stability: Conditions and Yes <input checked="" type="checkbox"/> No <input type="checkbox"/>	Materials to Avoid	Hazardous Decomposition Products: None.	
SECTION XII - TOXICOLOGICAL INFORMATION			
Route of Entry: <input checked="" type="checkbox"/> Skin Contact <input checked="" type="checkbox"/> Ingestion <input type="checkbox"/> Eye Contact <input type="checkbox"/> Inhalation Acute <input type="checkbox"/> Inhalation Chronic <input checked="" type="checkbox"/> Ingestion		Specific and Symptoms of Short-Term Exposure: None known.	
Chronic Effects and Symptoms of Long-Term Exposure: None known.		Sanitization: No sensitivity affects known.	
SECTION XIII - ECOLOGICAL INFORMATION		SECTION XIV - DISPOSAL CONSIDERATIONS	
Biodegradability: Not defined.		Method of Waste Disposal: Dispose in accordance with local governmental regulations.	
Further Details: N/A		SECTION XV - TRANSPORT INFORMATION	
Not regulated.		SECTION XVI - REGULATORY INFORMATION	
N/A.		SECTION XVII - OTHER INFORMATION	
FOR DENTAL USE ONLY: Use as directed. The information and recommendations are taken from sources (raw material MSDS's) and manufacturer's knowledge believed to be accurate; however, Ultradent Products, Inc., makes no warranty with respect to the accuracy of the information or the suitability of the recommendation and assumes no liability to any user thereof. Each user should review these recommendations in the specific context of the intended use and determine whether they are appropriate.			

© Copyright 2006 Ultradent Products, Inc. All Rights Reserved.

ULTRADENT
 PRODUCTS, INC.