Endo-Eze® Motor & Handpiece

Description & Indications:

The Endo-Eze® air motor (REF/UP 1178) is cleaned and packaged in a plastic box, and not sterilized in the delivery state. The Endo-Eze® Contra-Angle Handpiece (REF/UP 1155) is designed to be used with the Endo-Eze Handling Files and all other files with handle diameters of 2.35mm. Before initial use, we recommend sterilizing the air motor, the Contra-Angle and the rotary instrument.

Assembly:

1. Lubricate 1178

- Insert spray cap (with nozzle) onto spray can (Fig. 1). Remove blue insert in nozzle and replace with black cone shaped insert.
- Shake spray can several times.
- Set the speed ring of the air motor to maximum (position F or R).
- Fit the drive air tube (the smallest air intake) of motor onto the spray cap.
- Keep spray can upright (Fig. 2).
- Spray for approximately 1 second.

Keep spray cap with nozzle for use with subsequent oil cans. Following the application of oil and before first use of the air motor, sterilization in an autoclave is strongly recommended.

2. Lubricate 1155

- Screw E-type nozzle onto the spray cap over the arrow-head nozzle (Figs. 3 & 4). Turn approximately 10 times.
- Insert spray cap onto spray can.
- Shake spray can heavily several times.
- Insert the nozzle into the back end of contra-angle handpiece (Fig. 5).
- Hold the contra angle handpiece firmly.
- Keep spray can upright.
- Spray for approximately 1 second.

3. Endo Files

- The Contra-Angle handpiece has a push-button chuck. No tool is needed to change the root canal file.
- To insert file, press the button firmly with thumb and insert file into chuck until it is secure (Fig. 6). Make sure file is secure.
- To remove file, press button firmly with thumb and remove file.

4. Connect Air Motor 1178

- Align the motor back-end tubes with holes in the supply hose connector (Fig. 7).
- Push the motor into the supply hose connector.
- Turn the supply hose nut and securely fasten it. If loose, air and/or water leak may take place, and air pressure to motor will become insufficient.
- When motor speed feels low, check this connection first.

5. Assemble Contra-Angle Handpiece

- Firmly seat Contra-Angle onto air motor until it clicks in place (Fig. 8).

6. Test Run Handpiece

- Hold handpiece with head downward and let run approximately 30 seconds (without root canal file attached).
- Wipe off handpiece with gauze or cellulose paper and prepare for sterilization.

7. Test Run Air Motor

- Start the air motor.
- Screw speed regulation ring to "F" clockwise rotation or "R" counter-clockwise rotation.
- Vary the air motor from minimum to maximum within 30 seconds.
- Wipe off the air motor with gauze or cellulose paper.

8. To Remove Handpiece

- Remove handpiece from motor by pressing the push-button on the air motor (Fig. 9).

9. To Remove Air Motor

- Unscrew the supply hose nut.
- Pull out of the supply hose (Fig. 10)

WARNING! When removing the air motor, make sure that the speed regulating ring is set to position F or R to prevent it coming off due to any buildup of pressure. Do not remove the air motor during operation.

Operating condition and functional check

1. Check tightness

- Check tightness at the connection between air motor and supply hose. Leakage reduces the performance of the air motor and quality of spray.
- Check tightness of handpiece and air motor by running motor with spray on. If air and water come out from the side of the Contra-Angle, change O-rings on motor attachment.

2. Regulate water flow

- Adjust valves on the unit.

Disinfection/Sterilization

1. Manual disinfection

- Disinfect external surfaces by spraying or wiping with microbiologically tested surface disinfectant (Fig. 11).

2. Manual Cleaning

- Wear protective gloves.
- Make sure the speed regulating ring is set to position F or R.
- Remove root canal file from handpiece and remove handpiece from air motor. Disconnect the air motor from the supply hose (Fig. 12).
- Clean handpiece and air motor externally with cloth or brush soaked in disinfectant.

**If there is a high calcium content in the spray water (>18dH represents 180mg CaO in 1 quart, ~ 1 liter water), we recommend blowing out the water tubes before sterilization. For handpiece, blow through water tubes with air syringe.*

3. Sterilization

- We recommend autoclaving up to 135° C for 15 minutes.

The air motor should be cleaned and lubricated before placing into autoclave. Only sterilize with distilled water. For sterilization use sterilization packaging. We recommend sterilizing after every patient.

Oil Service

Keep away from open flames or ignition sources. Do not smoke near this product. Lubricate before every sterilization or twice daily after 30 minutes of continuous use.

Storage

Keep handpiece and air motor in sealed sterilization packaging (Fig. 13-14).

Description of the air motor

#1 Coupling, handpiece side with positive lock provided by push button 3.


#2 O-Rings. 3 O-rings Spray connection integrated.

#3 Push button lock

#4 Rotary direction. Clockwise and counter-clockwise direction selected with speed regulation ring:
F = Clockwise
R = Counter-clockwise

#5 Supply hose connection threads.

#6 Rubber gasket.

Operating pressure

- The operating pressure is measured as flow pressure.
- Operating pressure possible with 32psi (2.2 bar) to 43psi (3.0 bar).
- If only the air motor is used on the supply hose, we advise you to set the operating pressure to 43psi (3.0 bar).
- If the air motor is used alternately with a turbine on the same supply hose, it is essential to ensure that the pressure set for the two drives is within the prescribed range.

Changing the O-rings

Exchange damaged O-rings immediately!

- O-Rings can be exchanged without a special tool.
- Do not use a sharp tool.
- Firmly squeeze together the worn O-Ring between thumb and forefinger so as to form a loop and then withdraw the O-ring (Fig. 15).
- Slide on new O-rings by hand.

Fig. 15 Removing O-rings


Fig. 1 Spray cap with nozzle.


Fig. 2 Holding the air motor firmly, hold the oil can upright and spray for approximately 1 second.


Fig. 3 Oil spray can.

Fig. 4 Spray cap with E-type nozzle.


Fig. 5 Hold handpiece firmly. Hold spray can upright and spray for approximately 1 second.

Fig. 6 Insert and chuck the root canal file.


Fig. 7 Connect motor to supply hose. Securely tighten the supply hose nut.


Fig. 8 Assemble the contra-angle handpiece.


Fig. 9 Remove Contra-Angle handpiece by pressing push-button on motor.


Fig. 10 Remove supply hose.

Fig. 11 Disinfect with surface disinfectants or ethyl alcohol.


Fig. 12 To remove air motor.

Fig. 13 Handpiece in sterile goods packaging.


Fig. 14 Air motor in sterile goods packaging.

| Technical Data - Air Motor | 1178 | Technical Data - Contra-angle | 1155 |
|---|----------------------------------|--|-------------------------------|
| Rated speed (r.p.m.) at an 32psi (2.2 bar) operating pressure of: 43psi (3.0 bar) | 5000 - 20000 5000 - 25000 | Transmission Ratio | 4:1 |
| Stationary torque (Ncm) | 4 | Color coding | green ring |
| Max. power (W) | 20 | Motor coupling according to standard | ISO 3964 |
| Air consumption (NI/min) 32psi (2.2 bar) 43psi (3.0 bar) | 42 50 | Recommended root canal files | files with shank of ø 2.35 mm |
| Spray water volume (ml/min) according to ISO 13294 | >100 | Minimum chucking length | until fully located |
| Operating pressure psi (bar) | 32psi - 43psi (2.2 - 3.0bar) | Maximum rated reamer/file speed (min-1) | 6 000 |
| Exhaust air pressure psi (bar) | 7.15psi (< 0.5bar) | Spray water volume (mL/min) according to ISO 7785-2 | >50 |
| Water pressure psi (bar) | 21.8psi - 36.3psi (1.5 - 2.6bar) | Recommended water pressure (bar) | 28psi-43psi (2-3bar) |
| Spray air pressure psi (bar) must be higher than the water pressure | 28psi - 43psi (2 - 3bar) | Recommended chip air pressure psi (bar) must be higher than water pressure | 28psi-56psi (2-4bar) |

In case of service or repair, contact Ultradent Products, Inc.

U.S. federal law restricts this device to sale by or on the order of a dentist.

Keep out of reach of children.

For immediate reorder and/or complete descriptions of Ultradent's product line,


refer to Ultradent's catalog or call Toll Free 1-800-552-5512.

Outside U.S. call (801) 572-4200 or visit www.ultradent.com.


© Copyright 2003 Ultradent Products, Inc. All Rights Reserved.

Endo-Eze[®] Motor and Handpiece

DE-Motoreinheit und Winkelstück/FR-Moteur et pièce à main Endo-Eze/NL-motor en handstuk/IT-Motore e manipo/ES-Motor y Pieza de Mano/PT-Motor e Peça de Mão/SV-motor och insats/DA-motor og håndtag/FI-Moottori ja käsikappale/EL-Μικρομότορ και Χειρολαβή

English Lubricate prior to first use. Then twice daily or prior to sterilization. Motor: Set the speed ring to position F or R.

Assemble the contra-angle handpiece. To remove handpiece from motor press the push-button.

Insert file into head of handpiece. Make sure file is secure.

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


English Align the motor back-end tubes with holes in the supply hose connector. Push the motor into the supply hose connector. Turn the supply hose hub to secure.

Disinfect with surface disinfectants or ethyl alcohol.

Sterilization in autoclave up to 135° C and use sterilization packaging.

Precautions:
1. When removing the air motor, make sure that the line marking of the speed regulating ring is set to position F or R to prevent it slipping out from buildup of pressure.
2. Coupling is not sterilizable.
3. Operating pressure possible with 2.2 bar (32psi) to 3.0 bar (43psi).

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