

Astringedent®

Spot Remover

REF/UP 2160
30ml

UP® **ULTRADENT**
PRODUCTS, INC.

Contains: 25% phosphoric acid in aqueous solution.

Directions: Place spotted area of clothing on a hard surface. Express sufficient quantity of spot remover to soak entire area of spot. Scrub with soft toothbrush moistened with additional spot remover. Rinse thoroughly with lukewarm water.

Cautions: Contains acid. Handle with care. Wear latex gloves. Keep Astringedent® Spot Remover away from counter tops and skin. Do NOT use intraorally or internally. Keep out of reach of children. Do NOT mix with anything except water.



Hazard Rating

- 4 = Severe
- 3 = Serious
- 2 = Moderate
- 1 = Slight
- 0 = Minimal

MSDS Required

Manufactured by:
Ultradent Products, Inc.
505 West 10200 South
South Jordan, Utah 84095, USA
Toll Free 1.800.552.5512

www.ultradent.com

#64563.1 022107