



Safety Data Sheet

acc. to OSHA HCS

Printing date 08/15/2025

Reviewed on 08/15/2025

1 Identification

· **Product identifier**

· **Trade name:** FORMA

· **Article number:**

SDS 442-001.02R01, FO00002, FO00003, 4201-P2, 4201P2-P3, FO00004, FO00005, 4202-P2, 4202P2-P3, FO00006, FO00007, 4263-P2, 4264-P2, 4265-P2, FO00008, FO00009, 4265P2-P3, 4266-P2, FO00010, FO00011, 4266P2-P3, 4267-P2, FO00012, FO00013, 4267P2-P3, 4268-P2, FO00014, FO00015, 4268P2-P3, 4269-P2, FO00016, FO00017, 4269P2-P3, 4270-P2, FO00018, FO00019, 4270P2-P3, 4271-P2, FO00020, FO00021, 4271P2-P3, 4272-P2, FO00022, FO00023, 4272P2-P3, 4273-P2, FO00031, FO00032, 4273P2-P3, 4274-P2, FO00033, FO00034, 4274P2-P3, 4275-P2, FO00044, FO00037, 4275P2-P3, 4276-P2, FO00001, FO00038, 4276P2-P3, 4277-P2, FO00030, FO00040, 4277P2-P3, 4280-P2, FO00042, FO00026, 4281-P2, 42812-P2, FO00043, FO0000G, 42813-P2, 4281P2-P3, FO00046, FO00047, 4286-P2, 4286P2-P3, FO0000S, 4287-P2, 4287P2-P3, 4288-P2, 4288P2-P3, 4289-P2, 4289P2-P3, 4290-P2, 4290P2-P3, 4291-P2, 4291P2-P3, 4292-P2, 4292P2-P3, 4293-P2, 4293P2-P3, 4294-P2, 4294P2-P3, 4295-P2, 4295P2-P3, 4296-P2, 4296P2-P3, 4297-P2, 4297P2-P3, 4298-P2, 4299-P2, 4300-P2, 4301-P2, 4302-P2, 4303-P2, 4304-P2, 4304P2-P3, 4305-P2, 4305P2-P3, 4306-P2, 4306P2-P3, 4307-P2, 4307P2-P3

· **Application of the substance / the mixture** Professional Dental Composite

· **Details of the supplier of the safety data sheet**

· **Manufacturer/Supplier:**

Ultradent Products Inc.
505 W. Ultradent Drive (10200 S)
South Jordan, UT 84095-3942
USA
onlineordersupport@ultradent.com
(800) 552-5512

· **Information department:** Customer Service

· **Emergency telephone number:**

CHEMTREC (NORTH AMERICA) : +1 (800) 424-9300
(INTERNATIONAL) : +(703) 527-3887

2 Hazard(s) identification

· **Classification of the substance or mixture**



GHS08 Health hazard

Toxic to Reproduction 1B

H360 May damage fertility or the unborn child.

Specific Target Organ Toxicity - Repeated Exposure 2 H373 May cause damage to the lung through prolonged or repeated exposure. Route of exposure: Inhalation.



GHS07

Sensitization - Skin 1

H317 May cause an allergic skin reaction.

Eye Irritation 2B

H320 Causes eye irritation.

(Contd. on page 2)

Safety Data Sheet

acc. to OSHA HCS

Printing date 08/15/2025

Reviewed on 08/15/2025

Trade name: **FORMA**

(Contd. of page 1)

- **Label elements**

- **GHS label elements** Void
- **Hazard pictograms** GHS07, GHS08
- **Signal word** Danger

- **Health Hazard-determining components of labeling:**

Diurethane Dimethacrylate
Hydrophobic Pyrogenic Silica
Ethyl-4-Dimethylamino Benzoate
Triethylene Glycol Dimethacrylate

- **Hazard statements**

H320 Causes eye irritation.
H317 May cause an allergic skin reaction.
H360 May damage fertility or the unborn child.
H373 May cause damage to the lung through prolonged or repeated exposure. Route of exposure: Inhalation.

- **Precautionary statements**

P201 Obtain special instructions before use.
P202 Do not handle until all safety precautions have been read and understood.
P260 Do not breathe dust/fume/gas/mist/vapors/spray.
P264 Wash thoroughly after handling.
P272 Contaminated work clothing must not be allowed out of the workplace.
P280 Wear protective gloves/protective clothing/eye protection/face protection.
P302+P352 If on skin: Wash with plenty of water.
P305+P351+P338 If in eyes: Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue rinsing.
P308+P313 IF exposed or concerned: Get medical advice/attention.
P321 Specific treatment (see on this label).
P314 Get medical advice/attention if you feel unwell.
P333+P313 If skin irritation or rash occurs: Get medical advice/attention.
P337+P313 If eye irritation persists: Get medical advice/attention.
P363 Wash contaminated clothing before reuse.
P405 Store locked up.
P501 Dispose of contents/container in accordance with local/regional/national/international regulations.

- **Classification system:**

- **NFPA ratings (scale 0 - 4)**



- **HMIS-ratings (scale 0 - 4)**



3 Composition/information on ingredients

- **Chemical characterization: Mixtures**

- **Description:** Mixture of the substances listed below with nonhazardous additions.

(Contd. on page 3)

Safety Data Sheet

acc. to OSHA HCS

Printing date 08/15/2025

Reviewed on 08/15/2025

Trade name: **FORMA**

(Contd. of page 2)

· Dangerous components:		
7631-86-9	Silicon Dioxide, Chemically Prepared	≥5-<25%
72869-86-4	Diurethane Dimethacrylate	≥1-≤10%
109-16-0	Triethylene Glycol Dimethacrylate	≥1-≤10%
1314-23-4	Zirconium Dioxide	≥1-≤10%
68909-20-6	Hydrophobic Pyrogenic Silica	≥1-≤10%
10287-53-3	Ethyl-4-Dimethylamino Benzoate	≥0.1-<1%

4 First-aid measures

· Description of first aid measures

· General information:

Symptoms of poisoning may even occur after several hours; therefore medical observation for at least 48 hours after the accident.

· After inhalation:

This product is a thick paste, therefore inhalation is extremely unlikely.

Seek medical treatment in case of complaints.

Supply fresh air and to be sure call for a doctor.

In case of unconsciousness place patient stably in side position for transportation.

· After skin contact:

If skin irritation continues, consult a doctor.

Immediately wash with water and soap and rinse thoroughly.

· After eye contact:

Rinse opened eye for several minutes under running water. If symptoms persist, consult a doctor.

· After swallowing: If symptoms persist consult doctor.

· Information for doctor:

· **Most important symptoms and effects, both acute and delayed** No further relevant information available.

· **Indication of any immediate medical attention and special treatment needed**

No further relevant information available.

5 Fire-fighting measures

· Extinguishing media

· **Suitable extinguishing agents:** Use fire fighting measures that suit the environment.

· Special hazards arising from the substance or mixture

During heating or in case of fire poisonous gases are produced.

· Advice for firefighters

· Protective equipment:

Wear fully protective suit.

Mouth respiratory protective device.

6 Accidental release measures

· Personal precautions, protective equipment and emergency procedures

Wear protective equipment.

Mount respiratory protective device.

· **Environmental precautions:** Do not allow to enter sewers/ surface or ground water.

· Methods and material for containment and cleaning up:

Absorb with liquid-binding material (sand, diatomite, acid binders, universal binders, sawdust).

(Contd. on page 4)

Safety Data Sheet

acc. to OSHA HCS

Printing date 08/15/2025

Reviewed on 08/15/2025

Trade name: **FORMA**

(Contd. of page 3)

Dispose contaminated material as waste according to section 13.
 Ensure adequate ventilation.
 · **Reference to other sections**
 See Section 7 for information on safe handling.
 See Section 8 for information on personal protection equipment.
 See Section 13 for disposal information.

7 Handling and storage

- **Handling:**
- **Precautions for safe handling**
 Ensure good ventilation/exhaustion at the workplace.
 Open and handle receptacle with care.
 Prevent formation of aerosols.
- **Information about protection against explosions and fires:** Keep respiratory protective device available.
- **Conditions for safe storage, including any incompatibilities**
- **Storage:**
- **Requirements to be met by storerooms and receptacles:** No special requirements.
- **Information about storage in one common storage facility:** Not required.
- **Further information about storage conditions:**
 Store in a cool place.
 Protect from exposure to the light.
 See product labelling.
 Keep receptacle tightly sealed.
- **Specific end use(s)** Professional Dental Composite

8 Exposure controls/personal protection

- **Additional information about design of technical systems:** No further data; see section 7.
- **Control parameters**
- **Components with limit values that require monitoring at the workplace:**
 The following constituents are the only constituents of the product which have a PEL, TLV or other recommended exposure limit.
 At this time, the other constituents have no known exposure limits.

7631-86-9 Silicon Dioxide, Chemically Prepared

ACGIH Long-term value: 10* 5** mg/m³IDLH Short-term value: 3000 mg/m³

1314-23-4 Zirconium Dioxide

PEL Long-term value: 5 mg/m³
as ZrREL Short-term value: 10 mg/m³
Long-term value: 5 mg/m³
as ZrTLV Short-term value: 10 mg/m³
Long-term value: 5 mg/m³
as Zr; A4

- **Additional information:** The lists that were valid during the creation were used as basis.

(Contd. on page 5)

Safety Data Sheet

acc. to OSHA HCS

Printing date 08/15/2025

Reviewed on 08/15/2025

Trade name: **FORMA**

(Contd. of page 4)

- **Exposure controls**

- **Personal protective equipment:**

- **General protective and hygienic measures:**

Keep away from foodstuffs, beverages and feed.

Immediately remove all soiled and contaminated clothing.

Wash hands before breaks and at the end of work.

Store protective clothing separately.

Avoid contact with the eyes.

Avoid contact with the eyes and skin.

- **Breathing equipment:**

In case of brief exposure or low pollution use respiratory filter device. In case of intensive or longer exposure use respiratory protective device that is independent of circulating air.

- **Protection of hands:**



Protective gloves

The glove material has to be impermeable and resistant to the product/ the substance/ the preparation.

Due to missing tests no recommendation to the glove material can be given for the product/ the preparation/ the chemical mixture.

Selection of the glove material is based on consideration of the penetration times, rates of diffusion and the degradation

- **Material of gloves**

The selection of suitable gloves does not only depend on the material, but also on further marks of quality and varies from manufacturer to manufacturer. As the product is a preparation of several substances, the resistance of the glove material can not be calculated in advance and has therefore to be checked prior to the application.

- **Penetration time of glove material**

The exact breakthrough time has to be found out by the manufacturer of the protective gloves and has to be observed.

- **Eye protection:**



Tightly sealed goggles

- **Body protection:** Protective work clothing

9 Physical and chemical properties

- **Information on basic physical and chemical properties**

- **General Information**

- **Appearance:**

Form: Fluid

Color: According to product specification

- **Odor:** Characteristic

- **Odor threshold:** Not determined.

- **pH-value:** Not applicable (non-aqueous)

- **Change in condition**

Melting point/Melting range: Undetermined.

Boiling point/Boiling range: >250 °C

(Contd. on page 6)

Safety Data Sheet

acc. to OSHA HCS

Printing date 08/15/2025

Reviewed on 08/15/2025

Trade name: **FORMA**

(Contd. of page 5)

· Flash point:	>150 °C
· Flammability:	Not applicable.
· Auto igniting:	255 °C
· Decomposition temperature:	Not determined.
· Ignition temperature:	Product is not selfigniting.
· Danger of explosion:	Product does not present an explosion hazard.
· Explosion limits:	
Lower:	Not determined.
Upper:	Not determined.
· Vapor pressure:	Not determined.
· Density at 20 °C:	0.69538-1.83525 g/cm ³
· Bulk density:	695-1,835 kg/m ³
· Relative density	Not determined
· Vapor density	Not determined.
· Evaporation rate	Not determined.
· Solubility in / Miscibility with Water:	Not miscible or difficult to mix.
· Partition coefficient (n-octanol/water):	Not determined.
· Viscosity:	
Dynamic:	Not determined.
Kinematic:	Not determined
· Other information	No further relevant information available.

10 Stability and reactivity

- **Reactivity Stable**
- **Chemical stability**
- **Thermal decomposition / conditions to be avoided:** No decomposition if used according to specifications.
- **Possibility of hazardous reactions** No dangerous reactions known.
- **Conditions to avoid**
 - Light
 - Heat
- **Incompatible materials:** No further relevant information available.
- **Hazardous decomposition products:** No dangerous decomposition products known.

11 Toxicological information

- **Information on toxicological effects**
- **Acute toxicity:**

· LD/LC50 values that are relevant for classification:		
72869-86-4 Diurethane Dimethacrylate		
Oral	LD50	>5,000 mg/kg (rat)

(Contd. on page 7)

Safety Data Sheet

acc. to OSHA HCS

Printing date 08/15/2025

Reviewed on 08/15/2025

Trade name: **FORMA**

(Contd. of page 6)

109-16-0 Triethylene Glycol Dimethacrylate

Oral	LD50	>5,000 mg/kg (rat)
	LC50 Fish	16.4 mg/l (Fish) (Toxicity to fish)
Dermal	LD50	>2,000 mg/kg (mouse)

· **Primary irritant effect:**· **on the skin:** No irritant effect.· **on the eye:** Irritating effect.· **Sensitization:** Sensitization possible through skin contact.· **Additional toxicological information:**

The product shows the following dangers according to internally approved calculation methods for preparations:
Irritant

· **Carcinogenic categories**· **IARC (International Agency for Research on Cancer)**

13760-80-0	Ytterbium Trifluoride	3
128-37-0	Butylated Hydroxytoluene	3
1309-37-1	Ferric Oxide	3
14807-96-6	Talc (Mg ₃ H ₂ (SiO ₃) ₄)	3

· **NTP (National Toxicology Program)**

None of the ingredients is listed.

· **OSHA-Ca (Occupational Safety & Health Administration)**

None of the ingredients is listed.

12 Ecological information· **Toxicity**· **Aquatic toxicity:****72869-86-4 Diurethane Dimethacrylate**

EC50	>0.6 mg/kg (Algae)
Biodegradability	28 days (Aerobic) (Biodegradability testing)

109-16-0 Triethylene Glycol Dimethacrylate

EC50	>100 mg/kg (Algae)
Biodegradability	28 days (Aerobic) (Biodegradability testing)
Aqua toxicity	32 mg/l (daphnia) (No Observed Effect Concentration)

· **Persistence and degradability** No further relevant information available.· **Behavior in environmental systems:**· **Bioaccumulative potential** No further relevant information available.· **Mobility in soil** No further relevant information available.· **Additional ecological information:**· **General notes:**

Water hazard class 1 (Self-assessment): slightly hazardous for water

Do not allow undiluted product or large quantities of it to reach ground water, water course or sewage system.

· **Results of PBT and vPvB assessment**· **PBT:** Not applicable.· **vPvB:** Not applicable.

(Contd. on page 8)

Safety Data Sheet

acc. to OSHA HCS

Printing date 08/15/2025

Reviewed on 08/15/2025

Trade name: **FORMA**

(Contd. of page 7)

· **Other adverse effects** No further relevant information available.

13 Disposal considerations

· **Waste treatment methods**· **Recommendation:**

Do not allow product to reach sewage system.

Disposal should be in accordance with applicable regional, national and local laws and regulations.

Dispose of contents/container in accordance with international, federal, state, and local regulations.

· **Uncleaned packagings:**· **Recommendation:** Disposal must be made according to official regulations.

14 Transport information

· **UN-Number**· **DOT, IMDG, IATA**

Not Regulated

· **UN proper shipping name**· **DOT, IMDG, IATA**

Not Regulated

· **Transport hazard class(es)**· **DOT, ADN, IMDG, IATA**· **Class**

Not Regulated

· **Packing group**· **DOT, IMDG, IATA**

Not Regulated

· **Environmental hazards:**

Not Applicable.

· **Special precautions for user**

Not Applicable

· **Transport in bulk according to Annex II of
MARPOL 73/78 and the IBC Code**

Not Applicable.

· **UN "Model Regulation":**

Not Regulated

15 Regulatory information

· **Safety, health and environmental regulations/legislation specific for the substance or mixture**

No further relevant information available.

· **Sara**· **Section 355 (extremely hazardous substances):**

None of the ingredients is listed.

· **Section 313 (Specific toxic chemical listings):**

None of the ingredients is listed.

· **TSCA (Toxic Substances Control Act):**

7631-86-9

Silicon Dioxide, Chemically Prepared

ACTIVE

72869-86-4

Diurethane Dimethacrylate

ACTIVE

109-16-0

Triethylene Glycol Dimethacrylate

ACTIVE

1314-23-4

Zirconium Dioxide

ACTIVE

(Contd. on page 9)

Safety Data Sheet

acc. to OSHA HCS

Printing date 08/15/2025

Reviewed on 08/15/2025

Trade name: FORMA

(Contd. of page 8)

68909-20-6	Hydrophobic Pyrogenic Silica	ACTIVE
10287-53-3	Ethyl-4-Dimethylamino Benzoate	ACTIVE
· Hazardous Air Pollutants		
None of the ingredients is listed.		
· Proposition 65		
· Chemicals known to cause cancer:		
None of the ingredients is listed.		
· Chemicals known to cause reproductive toxicity for females:		
None of the ingredients is listed.		
· Chemicals known to cause reproductive toxicity for males:		
None of the ingredients is listed.		
· Chemicals known to cause developmental toxicity:		
None of the ingredients is listed.		
· Carcinogenic categories		
· EPA (Environmental Protection Agency)		
None of the ingredients is listed.		
· ACGIH Carcinogenicity (American Conference of Governmental Industrial Hygienists)		
1314-23-4	Zirconium Dioxide	A4
13760-80-0	Ytterbium Trifluoride	A4
128-37-0	Butylated Hydroxytoluene	A4
1309-37-1	Ferric Oxide	A4
14807-96-6	Talc (Mg3H2(SiO3)4)	A4
· NIOSH-Ca (National Institute for Occupational Safety and Health)		
None of the ingredients is listed.		
· Chemical safety assessment:		
Device is biocompatible when used as directed by dental professionals per ISO 10993-1		

16 Other information

This information is based on our present knowledge. However, this shall not constitute a guarantee for any specific product features and shall not establish a legally valid contractual relationship.

- **Department issuing SDS:** Environmental, Health, and Safety
- **Contact:** Customer Service
- **Date of preparation / last revision** 08/15/2025 / -
- **Abbreviations and acronyms:**
 - IMDG: International Maritime Code for Dangerous Goods
 - DOT: US Department of Transportation
 - IATA: International Air Transport Association
 - EINECS: European Inventory of Existing Commercial Chemical Substances
 - ELINCS: European List of Notified Chemical Substances
 - CAS: Chemical Abstracts Service (division of the American Chemical Society)
 - NFPA: National Fire Protection Association (USA)
 - HMIS: Hazardous Materials Identification System (USA)
 - LC50: Lethal concentration, 50 percent
 - LD50: Lethal dose, 50 percent
 - PBT: Persistent, Bioaccumulative and Toxic
 - vPvB: very Persistent and very Bioaccumulative
 - NIOSH: National Institute for Occupational Safety
 - OSHA: Occupational Safety & Health
 - TLV: Threshold Limit Value

(Contd. on page 10)

Safety Data Sheet
acc. to OSHA HCS

Printing date 08/15/2025

Reviewed on 08/15/2025

Trade name: FORMA

(Contd. of page 9)

*PEL: Permissible Exposure Limit**REL: Recommended Exposure Limit**Eye Irritation 2B: Serious eye damage/eye irritation – Category 2B**Sensitization - Skin 1: Skin sensitisation – Category 1**Toxic to Reproduction 1B: Reproductive toxicity – Category 1B**Specific Target Organ Toxicity - Repeated Exposure 2: Specific target organ toxicity (repeated exposure) – Category 2**. * Data compared to the previous version altered.*

US