



Safety Data Sheet acc. to OSHA HCS

Printing date 09/09/2026

Reviewed on 09/09/2026

1 Identification

- **Product identifier**
- **Trade name:** Opalescence™ Custom 35% Birthday Cake
- **Article number:** SDS 534-001.01R01, 7047-US, 1014097, 7179-US
- **Details of the supplier of the safety data sheet**
- **Manufacturer/Supplier:**
Ultradent Products Inc.
505 W. Ultradent Drive (10200 S)
South Jordan, UT 84095-3942
USA
onlineordersupport@ultradent.com
(800) 552-5512
- **Information department:** Customer Service
- **Emergency telephone number:**
CHEMTREC (NORTH AMERICA) : +1 (800) 424-9300
(INTERNATIONAL) : +(703) 527-3887

2 Hazard(s) identification

- **Classification of the substance or mixture**



GHS07

Skin Irritation 2 H315 Causes skin irritation.

Eye Irritation 2A H319 Causes serious eye irritation.

- **Additional information:** The full text of the classification and labeling elements mentioned is given in Section 16.

- **Label elements**

- **GHS label elements** The product is classified according to the Globally Harmonized System (GHS)

- **Hazard pictograms** GHS07

- **Signal word** Warning

- **Hazard statements**

H315 Causes skin irritation.

H319 Causes serious eye irritation.

- **Precautionary statements**

P264 Wash thoroughly after handling.

P280 Wear protective gloves / eye protection / face protection.

P302+P352 If on skin: Wash with plenty of water.

P321 Specific treatment (see on this label).

P305+P351+P338 If in eyes: Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue rinsing.

P332+P313 If skin irritation occurs: Get medical advice/attention.

P362+P364 Take off contaminated clothing and wash it before reuse.

P337+P313 If eye irritation persists: Get medical advice/attention.

(Contd. on page 2)

US

Safety Data Sheet

acc. to OSHA HCS

Printing date 09/09/2026

Reviewed on 09/09/2026

Trade name: Opalescence™ Custom 35% Birthday Cake

(Contd. of page 1)

· **Classification system:**· **NFPA ratings (scale 0 - 4)**

Health = 2

Fire = 0

Reactivity = 0

· **HMIS-ratings (scale 0 - 4)**

Health = *2

Fire = 0

Reactivity = 0

3 Composition/information on ingredients

· **Chemical characterization: Mixtures**· **Description:** Mixture of the substances listed below with nonhazardous additions.· **Dangerous components:**

56-81-5	Glycerin	≥10-<40%
124-43-6	Carbamide Peroxide	>10-≤25%
9003-01-4	Polyacrylic Acid	1-10%
7722-84-1	Hydrogen Peroxide	≥5-<8%
25322-68-3	Polyethylene Glycol	1-10%
1310-73-2	Sodium Hydroxide	≥1-<3%
7681-49-4	Sodium Fluoride	≥0.1-<1%

4 First-aid measures

· **Description of first aid measures**· **General information:** Immediately remove any clothing soiled by the product.· **After inhalation:** In case of unconsciousness place patient stably in side position for transportation.· **After skin contact:** Immediately wash with water and soap and rinse thoroughly.· **After eye contact:**

Rinse opened eye for several minutes under running water. If symptoms persist, consult a doctor.

· **After swallowing:** If symptoms persist consult doctor.· **Information for doctor:**· **Most important symptoms and effects, both acute and delayed** No further relevant information available.· **Indication of any immediate medical attention and special treatment needed**

No further relevant information available.

5 Fire-fighting measures

· **Extinguishing media**· **Suitable extinguishing agents:** Use fire fighting measures that suit the environment.· **Special hazards arising from the substance or mixture** No further relevant information available.· **Advice for firefighters**· **Protective equipment:** No special measures required.

US

(Contd. on page 3)

Safety Data Sheet

acc. to OSHA HCS

Printing date 09/09/2026

Reviewed on 09/09/2026

Trade name: Opalescence™ Custom 35% Birthday Cake

(Contd. of page 2)

6 Accidental release measures

- **Personal precautions, protective equipment and emergency procedures** Not required
- **Environmental precautions:** Do not allow to enter sewers/ surface or ground water.
- **Methods and material for containment and cleaning up:** Pick up mechanically.
- **Reference to other sections**
See Section 7 for information on safe handling.
See Section 8 for information on personal protection equipment.
See Section 13 for disposal information.

7 Handling and storage

- **Handling:**
- **Precautions for safe handling** See product labeling.
- **Information about protection against explosions and fires:** No special measures required.
- **Conditions for safe storage, including any incompatibilities**
- **Storage:**
- **Requirements to be met by storerooms and receptacles:** No special requirements.
- **Information about storage in one common storage facility:** Not required
- **Further information about storage conditions:** Keep receptacle tightly sealed.
- **Specific end use(s)** No further relevant information available.

8 Exposure controls/personal protection

- **Additional information about design of technical systems:** No further data; see section 7.
- **Control parameters**
- **Components with limit values that require monitoring at the workplace:**
The following constituents are the only constituents of the product which have a PEL, TLV or other recommended exposure limit.
At this time, the remaining constituent has no known exposure limits.

56-81-5 Glycerin

PEL	Long-term value: 15* 5** mg/m ³ mist; *total dust **respirable fraction
TLV	TLV withdrawn-insufficient data human occup. exp.

9003-01-4 Polyacrylic Acid

TWA	Short-term value: 0.05 mg/m ³
-----	--

7722-84-1 Hydrogen Peroxide

PEL	Long-term value: 1.4 mg/m ³ , 1 ppm
REL	Long-term value: 1.4 mg/m ³ , 1 ppm
TLV	Long-term value: 1.4 mg/m ³ , 1 ppm A3

25322-68-3 Polyethylene Glycol

WEEL	Long-term value: 10 mg/m ³ (H); MW>200
------	--

1310-73-2 Sodium Hydroxide

PEL	Long-term value: 2 mg/m ³
REL	Ceiling limit value: 2 mg/m ³

(Contd. on page 4)

Safety Data Sheet

acc. to OSHA HCS

Printing date 09/09/2026

Reviewed on 09/09/2026

Trade name: Opalescence™ Custom 35% Birthday Cake

(Contd. of page 3)

TLV	Ceiling limit value: 2 mg/m ³
7681-49-4 Sodium Fluoride	
PEL	Long-term value: 2.5 mg/m ³ as F
REL	Long-term value: 2.5 mg/m ³ as F
TLV	Long-term value: 2.5 mg/m ³ as F, A4; BEI

· **Ingredients with biological limit values:**

7681-49-4 Sodium Fluoride

BEI	2 mg/L Medium: urine Time: prior to shift Parameter: Fluoride (background, nonspecific)
	3 mg/L Medium: urine Time: end of shift Parameter: Fluoride (background, nonspecific)

· **Additional information:** The lists that were valid during the creation were used as basis.

· **Exposure controls**

· **Personal protective equipment:**

· **General protective and hygienic measures:**

Keep away from foodstuffs, beverages and feed.

Immediately remove all soiled and contaminated clothing.

Wash hands before breaks and at the end of work.

Avoid contact with the eyes and skin.

· **Breathing equipment:** Not required

· **Protection of hands:**



Protective gloves

The glove material has to be impermeable and resistant to the product/ the substance/ the preparation.

Due to missing tests no recommendation to the glove material can be given for the product/ the preparation/ the chemical mixture.

Selection of the glove material is based on consideration of the penetration times, rates of diffusion and the degradation

· **Material of gloves**

The selection of suitable gloves does not only depend on the material, but also on further marks of quality and varies from manufacturer to manufacturer. As the product is a preparation of several substances, the resistance of the glove material can not be calculated in advance and has therefore to be checked prior to the application.

· **Penetration time of glove material**

The exact breakthrough time has to be found out by the manufacturer of the protective gloves and has to be observed.

(Contd. on page 5)

US

Safety Data Sheet

acc. to OSHA HCS

Printing date 09/09/2026

Reviewed on 09/09/2026

Trade name: Opalescence™ Custom 35% Birthday Cake

(Contd. of page 4)

· **Eye protection:**

Tightly sealed goggles

· **Body protection:** Protective work clothing

9 Physical and chemical properties

· **Information on basic physical and chemical properties**· **General Information**· **Appearance:**

Form:	Gel
Color:	Clear
Odor:	Birthday Cake
Odor threshold:	Not determined

· **pH-value at 20 °C:** 5.4-7.4· **Change in condition**

Melting point/Melting range:	Undetermined
Boiling point/Boiling range:	Undetermined

· **Flash point:** Not applicable· **Flammability:** Not determined· **Decomposition temperature:** Not determined· **Ignition temperature:** Product is not selfigniting.· **Danger of explosion:** Product does not present an explosion hazard.· **Explosion limits:**

Lower:	1.5 Vol %
Upper:	15 Vol %

· **Vapor pressure:** Not applicable· **Vapor pressure at 50 °C:** 280 hPa· **Density:** Not determined· **Relative density** Not determined· **Vapor density** Not applicable· **Evaporation rate** Not applicable.· **Solubility in / Miscibility with**

Water: Soluble

· **Partition coefficient (n-octanol/water):** Not determined· **Viscosity:**

Dynamic:	Not applicable
Kinematic:	Not applicable

· **Other information** No further relevant information available.
 US
 (Contd. on page 6)

Safety Data Sheet

acc. to OSHA HCS

Printing date 09/09/2026

Reviewed on 09/09/2026

Trade name: Opalescence™ Custom 35% Birthday Cake

(Contd. of page 5)

10 Stability and reactivity

- **Reactivity** No further relevant information available.
- **Chemical stability**
- **Thermal decomposition / conditions to be avoided:** No decomposition if used according to specifications.
- **Possibility of hazardous reactions** No dangerous reactions known.
- **Conditions to avoid** No further relevant information available.
- **Incompatible materials:** No further relevant information available.
- **Hazardous decomposition products:** No dangerous decomposition products known.

11 Toxicological information

- **Information on toxicological effects**
- **Acute toxicity:**

· **LD/LC50 values that are relevant for classification:**

ATE (Acute Toxicity Estimate)

Oral	LD50	2,624-4,109 mg/kg
Dermal	LD50	27,959 mg/kg
Inhalative	LC50/4 h	200 mg/l

56-81-5 Glycerin

Oral	LD50	7,750 mg/kg (guinea pig) 4,100 mg/kg (mouse) 5,570 mg/kg (rat) 27,000 mg/kg (rabbit)
Dermal	LC50 Fish LD50	>5,000 mg/l (Fish) >21,900 mg/kg (rat) 10,000 mg/kg (rabbit)
Inhalative	LC50/4 h	>0.1425 mg/l (rat)

124-43-6 Carbamide Peroxide

Oral	LD50	>2,000 mg/kg (rat)
------	------	--------------------

9003-01-4 Polyacrylic Acid

Oral	LC50 Fish	580 mg/l (Fish)
------	-----------	-----------------

7722-84-1 Hydrogen Peroxide

Oral	LC50 Fish	16.4 mg/l (Fish)
------	-----------	------------------

25322-68-3 Polyethylene Glycol

Oral	LD50	19,600 mg/kg (guinea pig) 17,300 mg/kg (mouse) >10,000 mg/kg (rat)
Dermal	LC50 Fish LD50 LC50(Daphnia magna)	>100 mg/l (Fish) >20,000 mg/kg (rabbit) >10,000 mg/l (Water Flea) (Toxicity to aquatic invertebrates)

1310-73-2 Sodium Hydroxide

Oral	LD50 LC50 Fish	130-340 mg/kg (rat) 160 mg/l (Fish)
------	-------------------	--

(Contd. on page 7)

US

Safety Data Sheet

acc. to OSHA HCS

Printing date 09/09/2026

Reviewed on 09/09/2026

Trade name: Opalescence™ Custom 35% Birthday Cake

(Contd. of page 6)

Dermal	LD50	1,350 mg/kg (rabbit)
	Absolute lethal concentration	180 ppm (Fish)
7681-49-4 Sodium Fluoride		
Oral	LD50	52 mg/kg (mouse)
	LC50 Fish (static)	17 mg/l (Fish)
Dermal	LD50	175 mg/kg (rat)

- **Primary irritant effect:**

- **on the skin:** Irritant to skin and mucous membranes.

- **on the eye:** Irritating effect.

- **Sensitization:** No sensitizing effects known.

- **Additional toxicological information:**

The product shows the following dangers according to internally approved calculation methods for preparations:
Irritant

- **Carcinogenic categories**

· IARC (International Agency for Research on Cancer)		
9003-01-4	Polyacrylic Acid	3
7722-84-1	Hydrogen Peroxide	3
7681-49-4	Sodium Fluoride	3

- **NTP (National Toxicology Program)**

None of the ingredients is listed.

- **OSHA-Ca (Occupational Safety & Health Administration)**

None of the ingredients is listed.

12 Ecological information

- **Toxicity**

Aquatic toxicity:	
56-81-5 Glycerin	
EC50	>10,000 mg/kg (Bacteria)
9003-01-4 Polyacrylic Acid	
EC50	174 mg/kg (daphnia)
7722-84-1 Hydrogen Peroxide	
EC50	1.38 mg/l (Algae) 2.4 mg/l (daphnia)
1310-73-2 Sodium Hydroxide	
EC50	40.38 mg/kg (Water Flea)
7681-49-4 Sodium Fluoride	
EC50	272 mg/kg (Algae) 98 mg/kg (daphnia)
Algae Toxicity (static)	7 mg/l (Algae)

- **Persistence and degradability** No further relevant information available.

- **Behavior in environmental systems:**

- **Bioaccumulative potential** No further relevant information available.

- **Mobility in soil** No further relevant information available.

(Contd. on page 8)

Safety Data Sheet

acc. to OSHA HCS

Printing date 09/09/2026

Reviewed on 09/09/2026

Trade name: Opalescence™ Custom 35% Birthday Cake

(Contd. of page 7)

- **Additional ecological information:**
- **General notes:**
Water hazard class 2 (Self-assessment): hazardous for water
Do not allow product to reach ground water, water course or sewage system.
Danger to drinking water if even small quantities leak into the ground.
- **Results of PBT and vPvB assessment**
- **PBT:** Not applicable
- **vPvB:** Not applicable
- **Other adverse effects** No further relevant information available.

13 Disposal considerations

- **Waste treatment methods**
- **Recommendation:**
Dispose of contents/container in accordance with international, federal, state, and local regulations.
- **Uncleaned packagings:**
- **Recommendation:** Disposal must be made according to official regulations.
- **Recommended cleansing agent:** Water, if necessary with cleansing agents.

*

14 Transport information

- | | |
|---|-----------------|
| · UN-Number | |
| · DOT, IMDG, IATA | Not Regulated |
| · UN proper shipping name | |
| · DOT, IMDG, IATA | Not Regulated |
| · Transport hazard class(es) | |
| · DOT, ADN, IMDG, IATA | |
| · Class | Not Regulated |
| · Packing group | |
| · DOT, IMDG, IATA | Not Regulated |
| · Environmental hazards: | Not Applicable. |
| · Special precautions for user | Not Applicable |
| · Transport in bulk according to Annex II of MARPOL 73/78 and the IBC Code | Not Applicable. |
| · UN "Model Regulation": | Not Regulated |

15 Regulatory information

- **Safety, health and environmental regulations/legislation specific for the substance or mixture**
- **Sara**

· Section 355 (extremely hazardous substances):

7722-84-1 Hydrogen Peroxide

· Section 313 (Specific toxic chemical listings):

7757-79-1 Potassium Nitrate

(Contd. on page 9)

Safety Data Sheet

acc. to OSHA HCS

Printing date 09/09/2026

Reviewed on 09/09/2026

Trade name: Opalescence™ Custom 35% Birthday Cake

(Contd. of page 8)

· TSCA (Toxic Substances Control Act):

56-81-5	Glycerin	ACTIVE
124-43-6	Carbamide Peroxide	ACTIVE
9003-01-4	Polyacrylic Acid	ACTIVE
7722-84-1	Hydrogen Peroxide	ACTIVE
25322-68-3	Polyethylene Glycol	ACTIVE
1310-73-2	Sodium Hydroxide	ACTIVE
7681-49-4	Sodium Fluoride	ACTIVE

· Hazardous Air Pollutants

None of the ingredients is listed.

· Proposition 65
· Chemicals known to cause cancer:

None of the ingredients is listed.

· Chemicals known to cause reproductive toxicity for females:

None of the ingredients is listed.

· Chemicals known to cause reproductive toxicity for males:

None of the ingredients is listed.

· Chemicals known to cause developmental toxicity:

64-17-5 Ethyl Alcohol

· Carcinogenic categories
· EPA (Environmental Protection Agency)

None of the ingredients is listed.

· ACGIH Carcinogenicity (American Conference of Governmental Industrial Hygienists)

7722-84-1	Hydrogen Peroxide	A3
64-17-5	Ethyl Alcohol	A3
7681-49-4	Sodium Fluoride	A4

· NIOSH-Ca (National Institute for Occupational Safety and Health)

None of the ingredients is listed.

· Chemical safety assessment: A chemical safety assessment has not been carried out.

16 Other information

This information is based on our present knowledge. However, this shall not constitute a guarantee for any specific product features and shall not establish a legally valid contractual relationship.

· Department issuing SDS: Regulatory Compliance

· Contact: Customer Service

· Date of preparation / last revision 09/09/2026 / -

· Abbreviations and acronyms:

IMDG: International Maritime Code for Dangerous Goods

DOT: US Department of Transportation

IATA: International Air Transport Association

EINECS: European Inventory of Existing Commercial Chemical Substances

ELINCS: European List of Notified Chemical Substances

CAS: Chemical Abstracts Service (division of the American Chemical Society)

NFPA: National Fire Protection Association (USA)

HMIS: Hazardous Materials Identification System (USA)

LC50: Lethal concentration, 50 percent

LD50: Lethal dose, 50 percent

(Contd. on page 10)

Safety Data Sheet

acc. to OSHA HCS

Printing date 09/09/2026

Reviewed on 09/09/2026

Trade name: Opalescence™ Custom 35% Birthday Cake

(Contd. of page 9)

*PBT: Persistent, Bioaccumulative and Toxic**vPvB: very Persistent and very Bioaccumulative**NIOSH: National Institute for Occupational Safety**OSHA: Occupational Safety & Health**TLV: Threshold Limit Value**PEL: Permissible Exposure Limit**REL: Recommended Exposure Limit**BEI: Biological Exposure Limit**Skin Irritation 2: Skin corrosion/irritation – Category 2**Eye Irritation 2A: Serious eye damage/eye irritation – Category 2A**· * Data compared to the previous version altered.*

US