



## Safety Data Sheet acc. to OSHA HCS

Printing date 09/30/2026

Reviewed on 09/30/2026

### 1 Identification

· **Product identifier**

· **Trade name:** Opalescence™ Quick PF 45% - Opalescence™ Custom 45% - Finished Good

· **Article number:**

SDS 89-1000.02R01, 15073, 10817, 13691, 15292, 5346, 5359-AU, 5348, 5358-AU, 5358-U, 7094-U, 7095-U, 7108-U, 5359-CN, 5359-U, 5359-US, 6199-US, 6264-US, 6279-US, 6279-U, 6279-AU

· **Details of the supplier of the safety data sheet**

· **Manufacturer/Supplier:**

Ultradent Products Inc.  
505 W. Ultradent Drive (10200 S)  
South Jordan, UT 84095-3942  
USA  
onlineordersupport@ultradent.com  
(800) 552-5512

· **Information department:** Customer Service

· **Emergency telephone number:**

CHEMTREC (NORTH AMERICA) : +1 (800) 424-9300  
(INTERNATIONAL) : +(703) 527-3887

### 2 Hazard(s) identification

· **Classification of the substance or mixture**



GHS07

Skin Irritation 2      H315 Causes skin irritation.

Eye Irritation 2A      H319 Causes serious eye irritation.

Sensitization - Skin 1      H317 May cause an allergic skin reaction.

· **Additional information:** The full text of the classification and labeling elements mentioned is given in Section 16.

· **Label elements**

· **GHS label elements** The product is classified according to the Globally Harmonized System (GHS)

· **Hazard pictograms** GHS07

· **Signal word** Warning

· **Health Hazard-determining components of labeling:**

Oils, Peppermint

· **Hazard statements**

H315 Causes skin irritation.

H319 Causes serious eye irritation.

H317 May cause an allergic skin reaction.

· **Precautionary statements**

P261      Avoid breathing dust/fume/gas/mist/vapors/spray

P264      Wash thoroughly after handling.

P272      Contaminated work clothing must not be allowed out of the workplace.

P280      Wear protective gloves / eye protection / face protection.

(Contd. on page 2)

US

# Safety Data Sheet

acc. to OSHA HCS

Printing date 09/30/2026

Reviewed on 09/30/2026

**Trade name: Opalescence™ Quick PF 45% - Opalescence™ Custom 45% - Finished Good**

(Contd. of page 1)

- P302+P352 *If on skin: Wash with plenty of water.*  
 P305+P351+P338 *If in eyes: Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue rinsing.*  
 P362+P364 *Take off contaminated clothing and wash it before reuse.*  
 P333+P313 *If skin irritation or rash occurs: Get medical advice/attention.*  
 P321 *Specific treatment (see on this label).*  
 P337+P313 *If eye irritation persists: Get medical advice/attention.*  
 P363 *Wash contaminated clothing before reuse.*  
 P501 *Dispose of contents/container in accordance with local/regional/national/international regulations.*

· **Classification system:**

· **NFPA ratings (scale 0 - 4)**



· **HMIS-ratings (scale 0 - 4)**



## 3 Composition/information on ingredients

· **Chemical characterization: Mixtures**

· **Description:** Mixture of the substances listed below with nonhazardous additions.

· **Dangerous components:**

|            |                     |          |
|------------|---------------------|----------|
| 56-81-5    | Glycerin            | ≥5-≤10%  |
| 124-43-6   | Carbamide Peroxide  | ≥5-≤10%  |
| 7722-84-1  | Hydrogen Peroxide   | ≥1-<5%   |
| 9003-01-4  | Polyacrylic Acid    | 1-10%    |
| 25322-68-3 | Polyethylene Glycol | 1-10%    |
| 1310-73-2  | Sodium Hydroxide    | ≥1-<3%   |
| 7757-79-1  | Potassium Nitrate   | ≥1-<3%   |
| 7681-49-4  | Sodium Fluoride     | ≥0.1-<1% |
| 8006-90-4  | Oils, Peppermint    | ≥0.1-<1% |

## 4 First-aid measures

· **Description of first aid measures**

· **General information:** Immediately remove any clothing soiled by the product.

· **After inhalation:**

Supply fresh air and to be sure call for a doctor.

In case of unconsciousness place patient stably in side position for transportation.

· **After skin contact:** Immediately wash with water and soap and rinse thoroughly.

· **After eye contact:**

Rinse opened eye for several minutes under running water. If symptoms persist, consult a doctor.

· **After swallowing:** If symptoms persist consult doctor.

· **Information for doctor:**

· **Most important symptoms and effects, both acute and delayed** No further relevant information available.

(Contd. on page 3)

# Safety Data Sheet

acc. to OSHA HCS

Printing date 09/30/2026

Reviewed on 09/30/2026

**Trade name: Opalescence™ Quick PF 45% - Opalescence™ Custom 45% - Finished Good**

(Contd. of page 2)

- **Indication of any immediate medical attention and special treatment needed**  
No further relevant information available.

## 5 Fire-fighting measures

- **Extinguishing media**
- **Suitable extinguishing agents:** Use fire fighting measures that suit the environment.
- **Special hazards arising from the substance or mixture** No further relevant information available.
- **Advice for firefighters**
- **Protective equipment:** No special measures required.

## 6 Accidental release measures

- **Personal precautions, protective equipment and emergency procedures** Not required
- **Environmental precautions:**  
Dilute with plenty of water.  
Do not allow to enter sewers/ surface or ground water.
- **Methods and material for containment and cleaning up:**  
Absorb with liquid-binding material (sand, diatomite, acid binders, universal binders, sawdust).  
Dispose contaminated material as waste according to section 13.  
Ensure adequate ventilation.
- **Reference to other sections**  
See Section 7 for information on safe handling.  
See Section 8 for information on personal protection equipment.  
See Section 13 for disposal information.

## 7 Handling and storage

- **Handling:**
- **Precautions for safe handling**  
See product labeling.  
Ensure good ventilation/exhaustion at the workplace.  
Prevent formation of aerosols.
- **Information about protection against explosions and fires:** No special measures required.
- **Conditions for safe storage, including any incompatibilities**
- **Storage:**
- **Requirements to be met by storerooms and receptacles:** No special requirements.
- **Information about storage in one common storage facility:** Not required
- **Further information about storage conditions:** Keep receptacle tightly sealed.
- **Specific end use(s)** No further relevant information available.

## 8 Exposure controls/personal protection

- **Additional information about design of technical systems:** No further data; see section 7.
- **Control parameters**
- **Components with limit values that require monitoring at the workplace:**  
The following constituents are the only constituents of the product which have a PEL, TLV or other recommended exposure limit.  
At this time, the other constituents have no known exposure limits.

(Contd. on page 4)

# Safety Data Sheet

## acc. to OSHA HCS

Printing date 09/30/2026

Reviewed on 09/30/2026

**Trade name: Opalescence™ Quick PF 45% - Opalescence™ Custom 45% - Finished Good**

(Contd. of page 3)

**56-81-5 Glycerin**

|     |                                                                                       |
|-----|---------------------------------------------------------------------------------------|
| PEL | Long-term value: 15* 5** mg/m <sup>3</sup><br>mist; *total dust **respirable fraction |
| TLV | TLV withdrawn-insufficient data human occup. exp.                                     |

**7722-84-1 Hydrogen Peroxide**

|     |                                                      |
|-----|------------------------------------------------------|
| PEL | Long-term value: 1.4 mg/m <sup>3</sup> , 1 ppm       |
| REL | Long-term value: 1.4 mg/m <sup>3</sup> , 1 ppm       |
| TLV | Long-term value: 1.4 mg/m <sup>3</sup> , 1 ppm<br>A3 |

**9003-01-4 Polyacrylic Acid**

|     |                                          |
|-----|------------------------------------------|
| TWA | Short-term value: 0.05 mg/m <sup>3</sup> |
|-----|------------------------------------------|

**25322-68-3 Polyethylene Glycol**

|      |                                                      |
|------|------------------------------------------------------|
| WEEL | Long-term value: 10 mg/m <sup>3</sup><br>(H); MW>200 |
|------|------------------------------------------------------|

**1310-73-2 Sodium Hydroxide**

|     |                                          |
|-----|------------------------------------------|
| PEL | Long-term value: 2 mg/m <sup>3</sup>     |
| REL | Ceiling limit value: 2 mg/m <sup>3</sup> |
| TLV | Ceiling limit value: 2 mg/m <sup>3</sup> |

· **Additional information:** The lists that were valid during the creation were used as basis.

· **Exposure controls**

· **Personal protective equipment:**

· **General protective and hygienic measures:**

Keep away from foodstuffs, beverages and feed.

Immediately remove all soiled and contaminated clothing.

Wash hands before breaks and at the end of work.

Avoid contact with the eyes and skin.

· **Breathing equipment:**

In case of brief exposure or low pollution use respiratory filter device. In case of intensive or longer exposure use respiratory protective device that is independent of circulating air.

· **Protection of hands:**



Protective gloves

The glove material has to be impermeable and resistant to the product/ the substance/ the preparation.

Due to missing tests no recommendation to the glove material can be given for the product/ the preparation/ the chemical mixture.

Selection of the glove material is based on consideration of the penetration times, rates of diffusion and the degradation

· **Material of gloves**

The selection of suitable gloves does not only depend on the material, but also on further marks of quality and varies from manufacturer to manufacturer. As the product is a preparation of several substances, the resistance of the glove material can not be calculated in advance and has therefore to be checked prior to the application.

· **Penetration time of glove material**

The exact breakthrough time has to be found out by the manufacturer of the protective gloves and has to be observed.

(Contd. on page 5)

US

# Safety Data Sheet

acc. to OSHA HCS

Printing date 09/30/2026

Reviewed on 09/30/2026

Trade name: Opalescence™ Quick PF 45% - Opalescence™ Custom 45% - Finished Good

(Contd. of page 4)

## · Eye protection:



Tightly sealed goggles

## · Body protection: Protective work clothing

## 9 Physical and chemical properties

## · Information on basic physical and chemical properties

## · General Information

## · Appearance:

|                 |                |
|-----------------|----------------|
| Form:           | Gel            |
| Color:          | Colorless      |
| Odor:           | Mint           |
| Odor threshold: | Not determined |

· pH-value at 20 °C: 5.6-7.2

## · Change in condition

|                              |              |
|------------------------------|--------------|
| Melting point/Melting range: | Undetermined |
| Boiling point/Boiling range: | Undetermined |

· Flash point: Not applicable

· Flammability: Not applicable

· Decomposition temperature: Not determined

· Ignition temperature: Product is not selfigniting.

· Danger of explosion: Product does not present an explosion hazard.

## · Explosion limits:

|        |                |
|--------|----------------|
| Lower: | Not determined |
| Upper: | Not determined |

· Vapor pressure: Not determined

· Density at 20 °C: 1.27 g/cm<sup>3</sup>

· Relative density: Not determined

· Vapor density: Not determined

· Evaporation rate: Not determined

## · Solubility in / Miscibility with

Water: Fully miscible.

· Partition coefficient (n-octanol/water): Not determined

## · Viscosity:

|            |                |
|------------|----------------|
| Dynamic:   | Not determined |
| Kinematic: | Not determined |

· Other information: No further relevant information available.

## 10 Stability and reactivity

· Reactivity: No further relevant information available.

· Chemical stability

· Thermal decomposition / conditions to be avoided: No decomposition if used according to specifications.

(Contd. on page 6)

# Safety Data Sheet

acc. to OSHA HCS

Printing date 09/30/2026

Reviewed on 09/30/2026

Trade name: Opalescence™ Quick PF 45% - Opalescence™ Custom 45% - Finished Good

(Contd. of page 5)

- **Possibility of hazardous reactions** No dangerous reactions known.
- **Conditions to avoid** No further relevant information available.
- **Incompatible materials:** No further relevant information available.
- **Hazardous decomposition products:** No dangerous decomposition products known.

## 11 Toxicological information

- **Information on toxicological effects**
- **Acute toxicity:**

· **LD/LC50 values that are relevant for classification:**

### ATE (Acute Toxicity Estimate)

|            |          |                   |
|------------|----------|-------------------|
| Oral       | LD50     | 3,986-5,209 mg/kg |
| Dermal     | LD50     | 29,461 mg/kg      |
| Inhalative | LC50/4 h | 305 mg/l          |

### 56-81-5 Glycerin

|        |           |                                                                                               |
|--------|-----------|-----------------------------------------------------------------------------------------------|
| Oral   | LD50      | 7,750 mg/kg (guinea pig)<br>4,100 mg/kg (mouse)<br>5,570 mg/kg (rat)<br>27,000 mg/kg (rabbit) |
| Dermal | LC50 Fish | >5,000 mg/l (Fish)                                                                            |
|        | LD50      | >21,900 mg/kg (rat)<br>10,000 mg/kg (rabbit)                                                  |

### 124-43-6 Carbamide Peroxide

|      |      |                    |
|------|------|--------------------|
| Oral | LD50 | >2,000 mg/kg (rat) |
|------|------|--------------------|

### 7722-84-1 Hydrogen Peroxide

|      |           |                  |
|------|-----------|------------------|
| Oral | LC50 Fish | 16.4 mg/l (Fish) |
|------|-----------|------------------|

### 9003-01-4 Polyacrylic Acid

|      |           |                 |
|------|-----------|-----------------|
| Oral | LC50 Fish | 580 mg/l (Fish) |
|------|-----------|-----------------|

### 25322-68-3 Polyethylene Glycol

|        |                     |                                                                          |
|--------|---------------------|--------------------------------------------------------------------------|
| Oral   | LD50                | 19,600 mg/kg (guinea pig)<br>17,300 mg/kg (mouse)<br>>10,000 mg/kg (rat) |
| Dermal | LC50 Fish           | >100 mg/l (Fish)                                                         |
|        | LD50                | >20,000 mg/kg (rabbit)                                                   |
|        | LC50(Daphnia magna) | >10,000 mg/l (Water Flea) (Toxicity to aquatic invertebrates)            |

### 1310-73-2 Sodium Hydroxide

|        |                               |                      |
|--------|-------------------------------|----------------------|
| Oral   | LD50                          | 130-340 mg/kg (rat)  |
|        | LC50 Fish                     | 160 mg/l (Fish)      |
| Dermal | LD50                          | 1,350 mg/kg (rabbit) |
|        | Absolute lethal concentration | 180 ppm (Fish)       |

### 7757-79-1 Potassium Nitrate

|        |           |                                           |
|--------|-----------|-------------------------------------------|
| Oral   | LD50      | 3,015 mg/kg (rat)<br>1,901 mg/kg (rabbit) |
|        | LC50 Fish | 1,378 mg/l (Fish)                         |
| Dermal | LD50      | >5,000 mg/kg (rat)                        |

(Contd. on page 7)

# Safety Data Sheet

acc. to OSHA HCS

Printing date 09/30/2026

Reviewed on 09/30/2026

Trade name: Opalescence™ Quick PF 45% - Opalescence™ Custom 45% - Finished Good

(Contd. of page 6)

|                                   |                     |                     |
|-----------------------------------|---------------------|---------------------|
|                                   | LC50(Daphnia magna) | 490 mg/l (daphnia)  |
| <b>7681-49-4 Sodium Fluoride</b>  |                     |                     |
| Oral                              | LD50                | 52 mg/kg (mouse)    |
|                                   | LC50 Fish (static)  | 17 mg/l (Fish)      |
| Dermal                            | LD50                | 175 mg/kg (rat)     |
| <b>8006-90-4 Oils, Peppermint</b> |                     |                     |
| Oral                              | LD50                | 2,490 mg/kg (mouse) |
|                                   |                     | 2,426 mg/kg (rat)   |

· **Primary irritant effect:**

· **on the skin:** Irritant to skin and mucous membranes.

· **on the eye:** Irritating effect.

· **Sensitization:** Sensitization possible through skin contact.

· **Additional toxicological information:**

The product shows the following dangers according to internally approved calculation methods for preparations:  
Irritant

· **Carcinogenic categories**

· **IARC (International Agency for Research on Cancer)**

|           |                   |   |
|-----------|-------------------|---|
| 7722-84-1 | Hydrogen Peroxide | 3 |
| 9003-01-4 | Polyacrylic Acid  | 3 |
| 7681-49-4 | Sodium Fluoride   | 3 |

· **NTP (National Toxicology Program)**

None of the ingredients is listed.

· **OSHA-Ca (Occupational Safety & Health Administration)**

None of the ingredients is listed.

## 12 Ecological information

· **Toxicity**

· **Aquatic toxicity:**

**56-81-5 Glycerin**

|      |                          |
|------|--------------------------|
| EC50 | >10,000 mg/kg (Bacteria) |
|------|--------------------------|

**7722-84-1 Hydrogen Peroxide**

|      |                    |
|------|--------------------|
| EC50 | 1.38 mg/l (Algae)  |
|      | 2.4 mg/l (daphnia) |

**9003-01-4 Polyacrylic Acid**

|      |                     |
|------|---------------------|
| EC50 | 174 mg/kg (daphnia) |
|------|---------------------|

**1310-73-2 Sodium Hydroxide**

|      |                          |
|------|--------------------------|
| EC50 | 40.38 mg/kg (Water Flea) |
|------|--------------------------|

**7681-49-4 Sodium Fluoride**

|                         |                    |
|-------------------------|--------------------|
| EC50                    | 272 mg/kg (Algae)  |
|                         | 98 mg/kg (daphnia) |
| Algae Toxicity (static) | 7 mg/l (Algae)     |

· **Persistence and degradability** No further relevant information available.

· **Behavior in environmental systems:**

· **Bioaccumulative potential** No further relevant information available.

· **Mobility in soil** No further relevant information available.

(Contd. on page 8)

# Safety Data Sheet

acc. to OSHA HCS

Printing date 09/30/2026

Reviewed on 09/30/2026

Trade name: Opalescence™ Quick PF 45% - Opalescence™ Custom 45% - Finished Good

(Contd. of page 7)

- **Additional ecological information:**
- **General notes:**  
Water hazard class 1 (Self-assessment): slightly hazardous for water  
Do not allow undiluted product or large quantities of it to reach ground water, water course or sewage system.
- **Results of PBT and vPvB assessment**
- **PBT:** Not applicable
- **vPvB:** Not applicable
- **Other adverse effects** No further relevant information available.

## 13 Disposal considerations

- **Waste treatment methods**
- **Recommendation:**  
Dispose of contents/container in accordance with international, federal, state, and local regulations.
- **Uncleaned packagings:**
- **Recommendation:** Disposal must be made according to official regulations.
- **Recommended cleansing agent:** Water, if necessary with cleansing agents.

## 14 Transport information

- |                                                                                   |                 |
|-----------------------------------------------------------------------------------|-----------------|
| · <b>UN-Number</b>                                                                |                 |
| · <b>DOT, IMDG, IATA</b>                                                          | Not Regulated   |
| · <b>UN proper shipping name</b>                                                  |                 |
| · <b>DOT, ADN, IMDG, IATA</b>                                                     | Not Regulated   |
| · <b>Transport hazard class(es)</b>                                               |                 |
| · <b>DOT, ADN, IMDG, IATA</b>                                                     |                 |
| · <b>Class</b>                                                                    | Not Regulated   |
| · <b>Packing group</b>                                                            |                 |
| · <b>DOT, IMDG, IATA</b>                                                          | Not Regulated   |
| · <b>Environmental hazards:</b>                                                   | Not Applicable. |
| · <b>Special precautions for user</b>                                             | Not Applicable  |
| · <b>Transport in bulk according to Annex II of MARPOL 73/78 and the IBC Code</b> | Not Applicable. |
| · <b>UN "Model Regulation":</b>                                                   | Not Regulated   |

## 15 Regulatory information

- **Safety, health and environmental regulations/legislation specific for the substance or mixture**
- **Sara**

### · Section 355 (extremely hazardous substances):

|           |                   |
|-----------|-------------------|
| 7722-84-1 | Hydrogen Peroxide |
|-----------|-------------------|

### · Section 313 (Specific toxic chemical listings):

|           |                   |
|-----------|-------------------|
| 7757-79-1 | Potassium Nitrate |
|-----------|-------------------|

### · TSCA (Toxic Substances Control Act):

|          |                    |        |
|----------|--------------------|--------|
| 56-81-5  | Glycerin           | ACTIVE |
| 124-43-6 | Carbamide Peroxide | ACTIVE |

(Contd. on page 9)

# Safety Data Sheet

## acc. to OSHA HCS

Printing date 09/30/2026

Reviewed on 09/30/2026

**Trade name: Opalescence™ Quick PF 45% - Opalescence™ Custom 45% - Finished Good**

(Contd. of page 8)

|            |                     |        |
|------------|---------------------|--------|
| 7722-84-1  | Hydrogen Peroxide   | ACTIVE |
| 9003-01-4  | Polyacrylic Acid    | ACTIVE |
| 25322-68-3 | Polyethylene Glycol | ACTIVE |
| 1310-73-2  | Sodium Hydroxide    | ACTIVE |
| 7757-79-1  | Potassium Nitrate   | ACTIVE |
| 7681-49-4  | Sodium Fluoride     | ACTIVE |
| 8006-90-4  | Oils, Peppermint    | ACTIVE |

· **Hazardous Air Pollutants**

None of the ingredients is listed.

· **Proposition 65**

· **Chemicals known to cause cancer:**

None of the ingredients is listed.

· **Chemicals known to cause reproductive toxicity for females:**

None of the ingredients is listed.

· **Chemicals known to cause reproductive toxicity for males:**

None of the ingredients is listed.

· **Chemicals known to cause developmental toxicity:**

None of the ingredients is listed.

· **Carcinogenic categories**

· **EPA (Environmental Protection Agency)**

None of the ingredients is listed.

· **ACGIH Carcinogenicity (American Conference of Governmental Industrial Hygienists)**

|           |                   |    |
|-----------|-------------------|----|
| 7722-84-1 | Hydrogen Peroxide | A3 |
| 7681-49-4 | Sodium Fluoride   | A4 |

· **NIOSH-Ca (National Institute for Occupational Safety and Health)**

None of the ingredients is listed.

· **Chemical safety assessment:** A chemical safety assessment has not been carried out.

## 16 Other information

This information is based on our present knowledge. However, this shall not constitute a guarantee for any specific product features and shall not establish a legally valid contractual relationship.

· **Department issuing SDS:** Regulatory Compliance

· **Contact:** Customer Service

· **Date of preparation / last revision** 09/30/2026 / -

· **Abbreviations and acronyms:**

IMDG: International Maritime Code for Dangerous Goods

DOT: US Department of Transportation

IATA: International Air Transport Association

EINECS: European Inventory of Existing Commercial Chemical Substances

ELINCS: European List of Notified Chemical Substances

CAS: Chemical Abstracts Service (division of the American Chemical Society)

NFPA: National Fire Protection Association (USA)

HMIS: Hazardous Materials Identification System (USA)

LC50: Lethal concentration, 50 percent

LD50: Lethal dose, 50 percent

PBT: Persistent, Bioaccumulative and Toxic

vPvB: very Persistent and very Bioaccumulative

NIOSH: National Institute for Occupational Safety

OSHA: Occupational Safety & Health

TLV: Threshold Limit Value

PEL: Permissible Exposure Limit

(Contd. on page 10)

**Safety Data Sheet**  
**acc. to OSHA HCS**

Printing date 09/30/2026

Reviewed on 09/30/2026

**Trade name: Opalescence™ Quick PF 45% - Opalescence™ Custom 45% - Finished Good***REL: Recommended Exposure Limit**Skin Irritation 2: Skin corrosion/irritation – Category 2**Eye Irritation 2A: Serious eye damage/eye irritation – Category 2A**Sensitization - Skin 1: Skin sensitisation – Category 1**· \* Data compared to the previous version altered.*

(Contd. of page 9)

US